3M ESPE

Imprint™ 3 Quick Step Putty

Imprint™ 3 Putty Soft

VPS Impression Material

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3M ESPE AG
Dental Products
D-82229 Seefeld - Germany

3M ESPE
Dental Products
St. Paul, MN 55144-1000

3M, ESPE, Garant and Imprint are trademarks of 3M or 3M ESPE AG.
© 3M 2007. All rights reserved.

4400075363001

ENGLISH

Product Description

Imprint™ 3 impression materials, manufactured by or for 3M ESPE, are addition type silicones for one-step and two-step impression technique. Imprint™ 3 Quick Step Putty and Imprint™ 3 Putty Soft are kneadable materials for handmixing. All other materials are available in the Garant cartridge and are statically mixed and directly applied with the Garant™ Dispenser, manufactured for 3M ESPE.

Quick Setting Products:

Product Name	ISO 4823	ISO 4823 Consistency	Paste Color	Mixing Ratio (Vol.)	Mixing Device
Imprint™ 3 Quick Step Putty	Type 0	putty	green	1:1	manually
Imprint™ 3 Quick Step Light Body	Type 3	light-bodied	pink	1:1	Garant Dispenser
Imprint™ 3 Quick Step Regular Body	Type 3	light-bodied	pink	1:1	Garant Dispenser
Imprint™ 3 Quick Step Ultra-Regular Body	Type 2	medium-bodied	caramel	1:1	Garant Dispenser

Regular Setting Products:

Product Name	ISO 4823	ISO 4823 Consistency	Paste Color	Mixing Ratio (Vol.)	Mixing Device
Imprint™ 3 Putty Soft	Type 0	putty	burgundy	1:1	manually
Imprint™ 3 Light Body	Type 3	light-bodied	yellow	1:1	Garant Dispenser
Imprint™ 3 Regular Body	Type 3	light-bodied	yellow	1:1	Garant Dispenser
Imprint™ 3 Ultra-Regular Body	Type 2	medium-bodied	petrol	1:1	Garant Dispenser

These Instructions for Use should be kept for the duration of product use. For all additionally mentioned products please refer to the corresponding Instructions for Use.

Areas of Application

- All precision impressions (of e. g.: crown, bridge, inlay and onlay preparations)

When One-step Technique is used, Imprint 3 Quick Step Putty is especially suitable for impressions of one or two prepared units.

Precautionary Measures

3M ESPE MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

Preparation

- Block out undercuts or interdental spaces (as with any other impression material with increased shore hardness) sufficiently to aid impression removal from the teeth after setting, as otherwise impression removal from the mouth will be difficult or could lead to the extraction of natural teeth or prostheses.

Impression Tray:

All rigid impression trays generally used for precision impressions are suitable.

- For sufficient adhesion, apply VPS Tray Adhesive, manufactured by 3M ESPE, on the tray and allow it to dry completely, 5–15 minutes drying time are ideal.

Garant Dispenser:

- Insert the cartridge with the wash impression material into the Garant Dispenser.
- Prior to application, check that both Garant cartridge openings are not clogged.
- When using new Garant cartridges, extrude a small amount of paste without the mixing tip, until base paste and catalyst are evenly extruded.
- Attach a yellow Garant™ Mixing Tip, manufactured for 3M ESPE, to the cartridge.
- For intraoral application attach a yellow Garant Intraoral Tip onto the Garant Mixing Tip.
- If required: enlarge the exit opening for the intraoral tip by cutting it with a scalpel.
- Make sure that base paste and catalyst are mixed completely and are extruded in a uniform color.
 - Paste color: see "Product Description".
 - Material that has set inside the Garant Mixing Tip should not be extruded by force as this may damage the Garant cartridge and Garant Mixing Tip and result in the formation of leaks.

Retraction

Suitable retraction agents include aluminium hydroxide chloride, aluminium sulfate or ferric sulfate solutions.

- Keep the areas from which the impression is to be taken dry during impression.
- The use of retraction cords may be required for subgingival preparations.
- Prior to the impression, thoroughly remove residual retraction material by flushing and drying.

Dosage and Mixing

- Use the respective dosage spoon to take equal volumes of the base paste and catalyst paste Imprint 3 Putty from the container, whereby the dosage spoons must be filled level. The total volume depends on the impression tray size.
- Knead the pastes manually for about 30 seconds. The mixed paste must display an even color without striping.
 - Paste color: please see "Product Description".

Times

Product	Working Time at 23°C/74°F from start of mixing min:sec	Intraoral Syringe Time min:sec	Intraoral Setting Time min:sec
Quick Setting			
Imprint 3 Quick Step Putty	1:30	–	2:30
Imprint 3 Quick Step Light Body	1:30*	0:40	2:30
Imprint 3 Quick Step Regular Body			
Imprint 3 Quick Step Ultra-Regular Body			
Intraoral setting time for all combinations of quick setting materials and one-step technique:			2:30
Regular Setting			
Imprint 3 Putty Soft	2:00	–	3:30
Imprint 3 Light Body	2:00*	1:00	3:30
Imprint 3 Regular Body			
Imprint 3 Ultra-Regular Body			
Intraoral setting time for all combinations of regular setting materials and one-step technique:			3:30

* Start of mixing = Entry of paste into the mixing tip.

These are the processing times for products that have been kept at a temperature of 23°C/74°F and 20–80% relative humidity. These times are shortened due to higher temperatures of the product or prolonged if the temperature of the product is lower.

Impression Taking

One-step Technique:

Recommended product combinations	
Imprint 3 Quick Step Putty	Imprint 3 Quick Step Regular Body or Imprint 3 Quick Step Ultra-Regular Body
Imprint 3 Putty Soft	Imprint 3 Regular Body or Imprint 3 Ultra-Regular Body

- Fill Imprint 3 Quick Step Putty or Putty Soft into the impression tray prepared with adhesive.
- While or after the tray is being loaded, apply one of the higher viscous wash materials, Imprint 3 (Quick Step) Regular Body or Imprint 3 (Quick Step) Ultra-Regular Body around the cleaned and dried preparation from the bottom up. Keep the application nozzle permanently immersed in the paste to avoid the formation of air bubbles.
 - The lower viscous wash materials Imprint 3 Quick Step Light Body and Imprint 3 Light Body may be displaced.
 - Do not exceed tray material and syringe material processing times to avoid incomplete or inaccurate impressions.
- Slowly place the loaded tray in the mouth parallel to the long axes of the prepared teeth, and hold it in position passively.
- Remove the impression from the mouth after completed setting.

Two-step Technique:

Recommended product combinations	
Imprint 3 Putty Soft	Imprint 3 Light Body or Imprint 3 Regular Body
Imprint 3 Quick Step Putty	Imprint 3 Quick Step Light Body or Imprint 3 Quick Step Regular Body

- Fill the Putty material into the impression tray prepared with adhesive.
- Slowly insert the loaded tray into the mouth parallel to the long axes of the prepared teeth, and hold it in place without exerting pressure.
- Remove the impression from the mouth after completed setting, clean and dry.
- Excise undercuts and interdental septa etc. and cut spillways, if required.
- Clean and dry the impression thoroughly to provide for optimal coalescence of the two materials.
- Insert the wash impression material into the preliminary impression at the appropriate sites. If required, apply around the cleaned and dried preparation, while keeping the application nozzle permanently immersed in the paste to avoid the formation of bubbles.
- Replace the preliminary impression into the mouth and hold it in place.
- Remove the impression from the mouth after completed setting.

After Impression Taking

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution for the period of time recommended by the manufacturer. Do not leave the impression submerged in the disinfection solution for more than the specified amount of time as otherwise the impression may be damaged.
- After disinfection, clean the impression of any remaining disinfection solution by rinsing under running water for approx. 15 seconds.

Model Preparation

- The impression can be poured 30 minutes after removal from the mouth. There are no other time limitations for pouring.
- In order to attain a bubble-free model, rinse the impression with water shortly before pouring, remove excess water or use a silicone surfactant.
- Pour the impression with a standard dental die or model stone. The impressions can be silver or copper plated.

Cleaning

- Garant Dispenser: Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135°C/275°F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
- Impression Tray: The adhesive can be removed from re-usable impression trays using acetone.

Notes

- Completely remove the oxygen inhibition layer of composite materials, e. g. fillings or core buildups, it may impair or retard the setting of silicone impression materials.
- Keep the filled Garant Mixing Tip on top of the cartridge to serve as a closure until the next use. Removing and replacing the spent Garant Mixing Tip may lead to carry-over of paste and ensuing formation of clogs.
- Vinyl polysiloxane impression materials must not be combined with condensation type silicones or polyether impression materials. Even small traces can affect the setting process.
- Disposable latex gloves impair the setting of silicone impression materials. Vinyl gloves are more suitable.
- Under no circumstances may impressions come into contact with solvent-containing fluids. This could result in expansion and inaccurate stone models.

Storage and Stability

Store the product at 15–25°C/59–77°F.

Do not use after the expiration date.

Store set impressions at temperatures below 30°C/86°F.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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