Permadyne™ Garant™ 2:1

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Times	Mixing min:sec	Processing time from start of mixing* min:sec	Setting time from start of mixing** min:sec	Residence time in the mouth** min:sec
Permadyne Garant 2:1		2:00	5:30	3:30
Permadyne Penta H		2:30	6:00	3:30
Permadyne heavy-bodied consistency	0:30	2:30	5:30	3:00

The working times shown are in conformity with ISO 4823:2000. Higher temperatures will shorten the working time, while lower temperatures will prolong it. We recommend a working time of 01:45 min:sec and compliance with the indicated setting time from the start of mixing to assure an ideal impression result at all times under clinical conditions.

* Start of mixing = entry of paste into the mixing tip.
 ** For product combinations, longer times are valid.

Impression Taking

- ▶ If necessary for removal of the set impression from the mouth appropriately block-out deep undercuts in interdental areas.

One Step/Two Viscosities Technique

- ▶ Load the selected material onto the impression tray according to the corresponding manufacturer instructions.
- ▶ At the same time or subsequently, apply Permadyne Garant 2:1 around the dry preparation from bottom up, keeping the tip of the Intraoral Tip immersed in the material at all times, and in contact with the surface of the tooth.
If required, fill a separate Elastomer syringe directly from the mixing tip without attaching the Intraoral Tip.
- ▶ Without delay, place the loaded tray in the mouth and hold it in place without applying pressure until setting is complete.

To remedy initial adhesion ("setting the valve"), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior position. If this proves difficult it may be necessary to carefully blow some air or water between the impression and the gingiva.

Functional impression

- ▶ After taking a functional impression with a heavy-bodied or medium-bodied consistency polyether material, use Permadyne Garant 2:1 for precise impression taking (correction).

After Impression Taking

- ▶ Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene

- ▶ Place the impression in a standard disinfectant solution commonly used for impressions, e.g. Impresept™, manufactured for 3M ESPE, for the period of time recommended by the manufacturer, i.e., 10 min in the case of Impresept. Excessive disinfection may damage the impression.
- ▶ After disinfection, rinse the impression under running water for approx. 15 sec.
 * Impresept is not available in all countries.

Model Preparation

- ▶ Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression. To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring. **Do not use surfactants as these impair the quality of polyether impressions and are not necessary!**
- ▶ Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

Cleaning

- ▶ **Garant Dispenser:** Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135°C/275°F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
- ▶ **Impression Tray:** The adhesive can be removed from impression trays using acetone.

Notes

- **Turning or reattaching the mixing tip without subsequent impression taking may lead to carry-over of paste and clogging.**
- At storage temperatures below 18°C/64°F, the pastes may contract such way that precise dosing cannot be guaranteed even at room temperature.
- Direct exposure to sunlight and damp storage conditions may damage the impression.
- Polyether impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modelling.
- 3M ESPE polyether materials may only be combined with 3M ESPE polyethers. The working and setting times of the individual products must be strictly observed.

Precautionary Measures

If an adverse reaction is observed in a patient, discontinue use of the product in this patient.

3M ESPE MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

Storage and Stability

Store the product at 18-25°C/64-77°F. **Do not refrigerate!** Do not use after the expiration date. Store impressions dry and below 30°C/86°F in the dark.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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ENGLISH

Product Description

Permadyne™ Garant™ 2:1, manufactured by 3M ESPE, is a light-bodied consistency polyether impression material supplied in the Garant cartridge. The mixing ratio is 2 volumes base paste:1 volume catalyst.

☞ For details on the products mentioned below, please refer to the corresponding Instructions for Use. These Instructions for Use should not be discarded for the duration of product use.

Areas of Application

- Impression of inlay, onlay, crown, bridge, and veneer preparations
- Functional impressions

Preparation

Garant™ dispenser/cartridge, manufactured for 3M ESPE:

- ▶ Upon use of a new cartridge and with each new mixing process, start without attaching the mixing tip by extruding a small quantity of paste until both components emerge uniformly. Remove clogged paste, if necessary.
- ▶ Subsequently, install the Garant™ Mixing Tip White and Garant™ Intraoral Tip White, both manufactured for 3M ESPE.

Impression trays:

Any impression trays generally used for precision impressions are suitable.

- ▶ For sufficient adhesion, apply a thin layer of Polyether Adhesive, manufactured by 3M ESPE, to the tray and allow to dry completely (at least 30-60 sec, 15 min drying time is ideal).

Retraction

Solutions based on aluminium hydroxide chloride or aluminium sulphate, are suitable retraction agents.

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3M, ESPE, Garant, Impresept and Permadyne are trademarks of 3M or 3M ESPE AG.
© 3M 2006. All rights reserved.