Do not use after the expiration date.

ENGLISH
General Information:
Adper™ Scotchbond™ SE Self-etch adhesive, manufactured by 3M ESPE, is a simple to use, two-bottle, self-etch bonding adhesive with bottle A containing the aqueous primer and bottle B containing the adhesive. The primer contains a color indicator which indicates complete coverage of the preparation. The color disappears upon taking the adhesive, indicating activation of the acid and insuring that a proper etch is added to the substrate. Scotchbond SE is a two-step, self-etching adhesive that functions a wide range of applications. These include bonding to all classes of light curing direct restorations as well as procedures involving porcelain, composite, and metal and root surface desensitization.

Recommendations:
Adper Scotchbond SE Self-etch adhesive is cured by exposure to visible light. The light curing times instructed with this product assume the use of a 3M ESPE light curing unit, manufactured by 3M ESPE, or other dental visible curing light of comparable intensity. Curing lights should be checked often for proper output using a reliable metering system. Do not use for drying should be free of oil and water contaminants.

Indications:
Bonding between dentin/enamel and composite filling materials
Bonding between dentin/enamel and compomer filling materials
Intraoral repair of porcelain, composite and metal using light-cure composite resin
Desensitization of hypersensitive toothroot surface desensitization

Contraindications:
Adper Scotchbond SE Self-etch adhesive is suited for purely light-curing materials, but not for dual-curing or self-curing materials, and is not indicated for indirect applications including in veneer placement.

Precautions:
General information for Patients:
This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If any allergic reactions occur, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel:
This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic exposure, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts skin, remove and discard glove, wash hands immediately with soap and water and then thoroughly clean exposed skin with water. 3M ESPE MSDS can be obtained from www.3m.com/MSDS or contact your local subsidiary.

Preparation
Remove the cavity preparation slightly by spraying with water. Use 2-3 brief blows of air to dry the tooth. Use caution should you be dried to the extent that there is no water visible on tooth surfaces. **Over-drying may result in post-operative sensitivity.** Use a no-touch technique to tooth structure with a rubber cup and an air and fluoride-free cleaning paste (e.g. pumice paste), rinse thoroughly with water spray and dry.

Pulp Exposure:
Vitreond™ Plus Light Cure Glass Ionomer Liner/base, manufactured by 3M ESPE, is used to seal dentin and pulp exposure. Discard this brush. If pulp exposure has occurred, use a minimum amount of calcium hydroxide (e.g. Akinline™) before application of 3M ESPE Vitreond Plus Liner/base. Adper Scotchbond SE Self-etch adhesive will bond to Vitreond Plus Liner/base.

Instructions for bonding to dentin and cut enamel with a light cure adhesive:
1) Isolation: A rubber dam is the preferred method of isolation.

- Lightly dry tooth.
- Dispense 1 drop of Liquid A into one of the mix wells, and 1 drop of Liquid B into the second mix well.
- Wet brush tip with Liquid A. Apply to the entire bonding area so that a continuous red-colored layer is obtained on the surface. Discard this brush.
- Wet second brush tip with Liquid B, and scrub into the entire wetted surface of the bonding area. The red color will disappear quickly, indicating that the etching components have been activated. Continue scrubbing with moderate finger pressure for 20 seconds to ensure a proper etch.
- Air dry thoroughly for 10 seconds to evaporate water. Adhesive should remain in place and be shiny in appearance upon completion of this step.
- Re-coat brush with Liquid B, and apply second coat to the entire bonding surface. Lightly air thin adhesive layer to adjust film thickness/consistency.
- Light cure for 10 seconds.
- Restorative placement: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

Note: for bonding to uncatal enamel a separate etch using 3M ESPE Scotchbond™ etchant is recommended before application of Adper Scotchbond SE Self-etch adhesive.

- Instructions for bonding to composite:**
- Isolation: Rubber dam is the preferred method of isolation.
 - Roughen the existing material: Roughen the surface of existing composite using either a bur, diamond or a sandblast technique. Use of a coarse diamond or sandblast technique. Dispense 1 drop of Liquid A into one of the mix wells, and 1 drop of Liquid B into the second mix well.
 - Wet brush tip with Liquid A. Apply to the entire bonding area so that a continuous red-colored layer is obtained on the surface. Discard this brush.
 - Wet second brush tip with Liquid B, and scrub into the entire wetted surface of the bonding area. The red color will disappear quickly, indicating that the etching components have been activated. Continue scrubbing with moderate finger pressure for 20 seconds to ensure a proper etch.
 - Air dry thoroughly for 10 seconds to evaporate water. Adhesive should remain in place and be shiny in appearance upon completion of this step.
 - Re-coat brush with Liquid B, and apply second coat to the entire bonding surface. Lightly air thin adhesive layer to adjust film thickness/consistency.
 - Light cure for 10 seconds.
 - Restorative placement: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

- Instructions for metal or porcelain repair:**
- Isolation: Rubber dam is the preferred method of isolation.
 - Preparation: Clean the surface to be repaired with a slurry of plain flour of pumice. Rinse and dry thoroughly. Roughen the surface of existing metal or porcelain using either a bur, diamond or a sandblast technique. Be careful to remove all loose porcelain and leave the margin. Remove surface glaze 1mm beyond the margin.
 - Dispense 1 drop of Liquid A into one of the mix wells, and 1 drop of Liquid B into the second mix well.
 - Wet brush tip with Liquid A. Apply to the entire bonding area so that a continuous red-colored layer is obtained on the surface. Discard this brush.
 - Wet second brush tip with Liquid B, and scrub into the entire wetted surface of the bonding area. The red color will disappear quickly, indicating that the etching components have been activated. Continue scrubbing with moderate finger pressure for 20 seconds to ensure a proper etch.
 - Air dry thoroughly for 10 seconds to evaporate water. Adhesive should remain in place and be shiny in appearance upon completion of this step.
 - Re-coat brush with Liquid B, and apply second coat to the entire bonding surface. Lightly air thin adhesive layer to adjust film thickness/consistency.
 - Light cure for 10 seconds.
 - Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

- Instructions for root surface desensitization:**
- Lightly clean the root surface with flour of pumice. Rinse and blot dry.
 - Dispense 1 drop of Liquid A into one of the mix wells, and 1 drop of Liquid B into the second mix well.
 - Wet brush tip with Liquid A. Apply to the entire root surface so that a continuous red-colored layer is obtained on the surface. Discard this brush.
 - Wet second brush tip with Liquid B, and scrub into the entire wetted surface of the root surface area. The red color will disappear quickly, indicating that the etching components have been activated. Continue scrubbing with moderate finger pressure for 20 seconds to ensure a proper etch.
 - Air dry thoroughly for 10 seconds to evaporate water. Adhesive should remain in place and be shiny in appearance upon completion of this step.
 - Re-coat brush with Liquid B, and apply second coat to the entire bonding surface. Lightly air thin adhesive layer to adjust film thickness/consistency.
 - Light cure for 10 seconds.
 - Remove the oxygen inhibited layer with a moistened gauze.

Additional Notes:
It is extremely important that the adhesive Liquid B is scrubbed into the prepared surface to the extent of 20 second period. Scrubbing the adhesive only partially or shortening the application will compromise the etch pattern and have a negative effect on adhesion. Non-free products should be used as, resins materials that contain eugenol, can adversely affect the polymerization of the adhesive.

- Storage and Shelf Life**
Adper Scotchbond SE Self-etch adhesive should be protected from light. Product should be stored between 2 °C/35 °F and 27°/C/80° F. Shelf life of Liquids A and B is 24 months at room temperature.

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