

RelyX™ Luting Plus Cement




RelyX™
Cements

RelyX™ Luting Plus Cement

- Convenient paste-paste delivery in the Clicker™ Dispenser
- Consistent mix ratios
- Enhanced formula for better performance
- Easy clean-up
- Low post-operative sensitivity
- Sustained fluoride release

3M ESPE, the leaders with a history of excellence in dental cementation, bring you the newest member of the RelyX cement family. RelyX™ Luting Plus Cement is designed to make your job easier with our convenient, paste-paste Clicker dispenser. The advanced resin-modified glass ionomer formula is easy to mix, load and seat, offers enhanced marginal integrity with low solubility, high fracture toughness, sustained fluoride release and virtually no post-operative sensitivity.

DATA: Shear Bond Strength Comparison



RelyX™ Luting Plus Cement provides increased adhesion and shear bond strength on dentin and enamel cementations.

RelyX™ Cements. Easy to choose. Easy to use.



RelyX™ Luting Plus Cement—the paste-paste formula performs even better than powder-liquid cements.



What Dentists are Saying

More and more dentists are discovering the convenience and performance of RelyX™ Luting Plus Cement. They applaud the Clicker™ Dispenser and paste-paste formula from 3M ESPE.

“We really like the ease of dispensing, mixing and especially the easy clean-up after seating.”


“We always loved RelyX™ Luting Cement and this delivery system just makes it better.”

“It was easy to dispense the right ratio every time. Clean-up was easier and the ability to release fluoride made this product even more appealing.”

Indications for Use

- Metal crowns, inlays and onlays (includes stainless steel crowns)
- Porcelain-fused-to-metal crowns and bridges
- All-zirconia or all-alumina strengthened core ceramic restorations (such as Lava™ or Procera®)
- Endodontic pins and posts
- Orthodontic appliances

Physical Property Comparison



Ordering Information

Item Number

3525TK



Product Information

RelyX™ Luting Plus Cement Trial Kit

1 - 11g RelyX™ Luting Plus Cement in Clicker™ Dispenser, illustrated technique guide and instructions for use.

3525



RelyX™ Luting Plus Cement Refills

2 - 11g RelyX™ Luting Plus Cement in Clicker™ Dispensers and instructions for use.

A Complete Family of Solutions

RelyX™ Cements are a complete family of high-quality, reliable products for all your cementation needs. For more information on all of the RelyX cements from 3M ESPE, visit our Web site, or call our technical hotline:

3M ESPE Technical Hotline: 1-800-634-2249

www.3MESPE.com


3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
USA

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-265-1840 ext. 6229



Printed in USA
© 3M 2005

3M, ESPE and RelyX and Lava are trademarks of 3M or 3M ESPE AG. FujiCEM and Fuji I are trademarks of GC America Inc. Procera is a trademark of Nobel Biocare. Used under license in Canada.

70-2009-3614-7